


Report Test

M/C : Henan Sanqiang Medical Equipment Liability Co. Ltd
Model : SQ-H220
Serial No : HNSQ2020050063

LOT No : 00963
Date : 08 April 2025

Determination of Sterilization Operation Specification

Item Detail	Specification	Setting	Result
Chamber Setting Temperature	40.0-55.5 °C	45.0 °C	
Vacuum pressure	40-60 kPa	60 kPa	
Vacuum Time	50-60 min	60 min	
Humidification			
Humidity	50.0-90.0%	20.2 %	
Sterilization			
Sterilization Temperature	40.0-55.5 °C	43.5 °C	
Gas Injection			
Final pressure	0-10 kPa	4 kPa	
Final temperature of vaporizer	40.0-55.5 °C	44.8 °C	
Exposure			
Temperature at start point of Exposure time period	40.0-55.5 °C	43.3 °C	
Temperature at finish point of Exposure time period	40.0-55.5 °C	46.0 °C	
Exposure Time	120-240 min	180 min	
Exhaust			
Exhaust Time	60-80 min	75 min	
First Vacuum pressure	0-2 kPa	1.5 kPa	
Second Vacuum pressure	0-2 kPa	1.5 kPa	
Third Vacuum pressure	0-2 kPa	1.5 kPa	
Fourth Vacuum pressure	0-2 kPa	1.5 kPa	
Fifth Vacuum pressure	0-2 kPa	1.5 kPa	
Storage			
Storage Time	≥ 48 hour	≥ 48 hour	
Chemical Test			
EO sterilization temperature 40-55°C, Time Sterilization by EO 3hours	Blue	Red	Pass
Biological Test			
EZ Test Biological Indicator	Green	Green	Pass


 (.....)
 Mr. Jutha Liammukda
 Research and Development Engineer



FM-ISO0-009 Rev.02 (05/11/2020)

			LOT	Color	BOX	TOTAL / PCS
		207379				
2	Z-GPUNO-1.3-6-2-35-CR	STERILE GOLD PLATED SURGICAL STEEL JEWELLED CURVED NOSE STUD	LOT NO 00963			100
3	Z-GPVLBJ-1.2-8-0.3-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWELLED MICRO LABRET	LOT NO 00963			100
		207404				
1	Z-GPVLBJ-1.2-10-0.3-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWELLED MICRO LABRET	LOT NO 00963			100
		207963				
23	Z-RGSCS413-CR	STERILE ROSE GOLD PVD COATED SURGICAL STEEL JEWELLED EAR STUDS PAIR - FLOWER	LOT NO 00963			5
		208010				
3	Z-SCBYC03-1.6	STERILE SURGICAL STEEL JEWELLED BELLY CLIPPER - FEATHER	LOT NO 00963			2
		208061				
1	Z-TNM/LBJD-1.2-10-0.3-HP-CR	STERILE TITANIUM INTERNALLY THREADED JEWELLED DISC MICRO LABRET	LOT NO 00963			15
2	Z-TNM/LBJD-1.2-8-0.3-HP-CR	STERILE TITANIUM INTERNALLY THREADED JEWELLED DISC MICRO LABRET	LOT NO 00963			15
		208062				
128	Z-BKTMN-1.2-8-0.3	STERILE BLACK PVD COATED TITANIUM CURVED MICRO BARBELL	LOT NO 00963			30
129	Z-BKTMN-1.2-8-0.3	STERILE BLACK PVD COATED TITANIUM MICRO LABRET	LOT NO 00963			30
133	Z-GPVLBJD-1.2-8-0.3-35-CR	STERILE GOLD PVD COATED TITANIUM PREMIUM CRYSTAL JEWELLED CURVED NOSE STUD	LOT NO 00963			5
		208067				
3	Z-BL-1.6-16-0.4	STERILE SURGICAL STEEL BARBELL	LOT NO 00963			40
4	Z-BL-1.6-16-0.5	STERILE SURGICAL STEEL BARBELL	LOT NO 00963			15
5	Z-BL-1.6-38-0.5	STERILE SURGICAL STEEL BARBELL	LOT NO 00963			5
5	Z-BN-1.6-10-0.5-8-CR	STERILE SURGICAL STEEL DOUBLE JEWELLED NAVEL BANANA	LOT NO 00963			10
7	Z-BN-1.6-10-0.5-8	STERILE SURGICAL STEEL NAVEL BANANA	LOT NO 00963			10
3	Z-BN-1.3-6-2-35-CR	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWELLED CURVED NOSE STUD	LOT NO 00963			30
9	Z-BN-1.2-10-0.4	STERILE SURGICAL STEEL CURVED MICRO BARBELL	LOT NO 00963			10
10	Z-SES1-6-CR	STERILE SURGICAL STEEL ROUND PROX-Q SET JEWELLED EAR STUDS PAIR	LOT NO 00963			20
		208068				
1	Z-MBN-1.2-10-0.4-CR	STERILE SURGICAL STEEL DOUBLE JEWELLED CURVED MICRO BARBELL	LOT NO 00963			20
2	Z-PFH-1.2-8-0.3	STERILE SURGICAL STEEL THREADLESS LABRET	LOT NO 00963			30
		208081				
1	Z-TMB-1.2-14-0.3-HP	STERILE TITANIUM MICRO BARBELL	LOT NO 00963			30
		208140				
1	Z-MCB-1.2-10-0.3	STERILE SURGICAL STEEL MICRO CIRCULAR BARBELL	LOT NO 00963			25
2	Z-COC-1.2-10-0.3	STERILE SURGICAL STEEL SEAMLESS RING	LOT NO 00963			10
		2532668				
1	Z-GPMLSSC377-1.2-8-0.3-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWELLED TRAGUS MICRO BARBELL - SNOWFLAKE	LOT NO 00963			20
		2534030				
14	PT-Z-TNMBN-1.2-8-0.3-HP	STERILE TITANIUM INTERNALLY THREADED CURVED MICRO BARBELL	LOT NO 00963			100
		2535417				
17	THD-Z-INSJ278-1.2-AM	STERILE SURGICAL STEEL JEWELLED MICRO ATTACHMENT FOR 1.2MM INTERNALLY THREADED PINS	LOT NO 00963			280
		2539450				
3	QL-Z-GPBN-1.6-12-0.5-8-CR	STERILE GOLD PVD COATED SURGICAL STEEL DOUBLE PREMIUM CRYSTALS JEWELLED NAVEL BANANA FOR QUILLAN	LOT NO 00963			50
		2539450				
2	TDI-Z-10G18TMLB-1.2-10-0.3	STERILE GOLD PVD 18K COATED TITANIUM MICRO LABRET	LOT NO 00963	SET/25		250
		2537984				
1	Z-TNMBN-1.2-12-0.3-HP	STERILE TITANIUM INTERNALLY THREADED CURVED MICRO BARBELL	LOT NO 00963			10
		2539080				
4	TDI-Z-10G18TML-1.6-14-0.5	STERILE GOLD PVD 18K COATED TITANIUM BARBELL	LOT NO 00963	SET/10		100
5	TDI-Z-10G18TML-1.6-15-0.5	STERILE GOLD PVD 18K COATED TITANIUM BARBELL	LOT NO 00963	SET/10		100
10	TDI-Z-10G18TMLB-1.2-8-0.3	STERILE GOLD PVD 18K COATED TITANIUM MICRO LABRET	LOT NO 00963	SET/30		300
		2538332				
1	PT-Z13-TMBN-1.2-8-0.3	STERILE TITANIUM MICRO CURVED BARBELL FOR PTT ROC	LOT NO 00963	SET/25		250
		2538406				
9	PT-Z-TBL-1.6-36-0.5	STERILIZED TITANIUM BARBELL FOR PTT ROC	LOT NO 00963			50
15	PT-Z13-TBN-1.6-10-0.5-8-CR	STERILE TITANIUM DOUBLE JEWELLED NAVEL BANANA FOR PTT ROC	LOT NO 00963	SET/10		400
16	PT-Z13-TBN-1.6-10-0.5-8-CR	STERILE TITANIUM JEWELLED NAVEL BANANA FOR PTT ROC	LOT NO 00963	SET/30		300
17	PT-Z13-TMBN-1.2-10-0.3	STERILE TITANIUM MICRO CURVED BARBELL FOR PTT ROC	LOT NO 00963	SET/25		250
10	PT-Z13-TMLB-1.2-8-0.3	STERILE TITANIUM MICRO LABRET FOR PTT ROC	LOT NO 00963	SET/100		1000
10	PT-Z13-TMLB-1.2-8-0.3	STERILE TITANIUM MICRO LABRET WITH JEWELLED BALL FOR PTT ROC	LOT NO 00963	SET/15		150
		2538468				
1	PT-Z13-TMBN-1.2-8-0.3	STERILE TITANIUM MICRO CURVED BARBELL FOR PTT ROC	LOT NO 00963	SET/25		250
		2538582				
3	TDI-Z-TMBN-1.2-10-0.3	STERILE TITANIUM MICRO CURVED BARBELL FOR TDI	LOT NO 00963			300
4	TDI-Z-TMLB-1.2-10-0.3	STERILE TITANIUM MICRO LABRET FOR TDI	LOT NO 00963			300
5	TDI-Z-BL-1.6-11-0.5	STERILE SURGICAL STEEL BARBELL FOR TDI	LOT NO 00963			200
		2537444				
5	SV-STERILIZE	SERVICE STERILE	LOT NO 00963			10
						5327

Sterilization:
Temperature: 44.3°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 15:12

Add FO gas:
Temperature: 44.3°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 15:12

Evacuation:
Temperature: 44.4°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 15:13

Keeping Temp:
Temperature: 45.2°C
Humidity: 20.38%
Pressure: -1.00Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 15:50

Keep pressure qualified:
Temperature: 43.7°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 15:54

Evacuation:
Temperature: 42.2°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:03

Heating-up:
Temperature: 20.4°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:20

Setting:
Temperature: 45.0°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:25

8/4/25
10 10 16.5

Air replacement:
Temperature: 45.3°C
Humidity: 20.38%
Pressure: -1.00Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:30

Air replacement:
Temperature: 45.3°C
Humidity: 20.38%
Pressure: -1.00Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:34

Air replacement:
Temperature: 45.4°C
Humidity: 20.38%
Pressure: -1.00Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:39

Air replacement:
Temperature: 45.1°C
Humidity: 20.38%
Pressure: -1.00Pa
Sterilization time: 3hour
Air replace rounds: 4
round time: 15min(s)
2025-4-8 16:39

Air replacement:
Temperature: 45.5°C
Humidity: 20.38%
Pressure: -1.00Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:41

Vacuumation:
Temperature: 45.5°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:41

Sterilization:
Temperature: 44.5°C
Humidity: 20.38%
Pressure: -0.90Pa
Sterilization time: 3hour
Air replace rounds: 5
round time: 15min(s)
2025-4-8 16:44

