

Report Test

M/C : Henan Sanqiang Medical Equipment Liability Co. Ltd
 Model : SQ-H220
 Serial No : HNSQ2020050063

LOT No : 00972
 Date : 12 May 2025

Determination of Sterilization Operation Specification

Item Detail	Specification	Setting	Result
Chamber Setting Temperature	40.0-55.5 °C	45.0 °C	
Vacuum pressure	40-60 kPa	60 kPa	
Vacuum Time	50-60 min	60 min	
Humidification			
Humidity	50.0-90.0%	20.2 %	
Sterilization			
Sterilization Temperature	40.0-55.5 °C	43.5 °C	
Gas Injection			
Final pressure	0-10 kPa	4 kPa	
Final temperature of vaporizer	40.0-55.5 °C	44.8 °C	
Exposure			
Temperature at start point of Exposure time period	40.0-55.5 °C	43.3 °C	
Temperature at finish point of Exposure time period	40.0-55.5 °C	46.0 °C	
Exposure Time	120-240 min	180 min	
Exhaust			
Exhaust Time	60-80 min	75 min	
First Vacuum pressure	0-2 kPa	1.5 kPa	
Second Vacuum pressure	0-2 kPa	1.5 kPa	
Third Vacuum pressure	0-2 kPa	1.5 kPa	
Fourth Vacuum pressure	0-2 kPa	1.5 kPa	
Fifth Vacuum pressure	0-2 kPa	1.5 kPa	
Storage			
Storage Time	≥ 48 hour	≥ 48 hour	
Chemical Test			
EO sterilization temperature 40-55°C, Time Sterilization by EO 3hours	Blue	Red	Pass
Biological Test			
EZ Test Biological Indicator	Green	Green	Pass

(*Jutha Liam*)
 Mr.Jutha Liamukda
 Research and Development Engineer

			LOT	Color	BOX	TOTAL / PCS
1	Z-XINLB-1.2-8.0-3-CL	STERILE BIOFLEX INTERNAL LABRET WITH PUSH FIT BALL	208587	LOT NO 00972		25
1	Z-GPINMLBJ109-1.2-8.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED JEWEL MICRO LABRET	208684	LOT NO 00972		6
2	Z-GPINMLBJ146-1.2-8.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED JEWEL MICRO LABRET		LOT NO 00972		4
3	Z-GPINMLBJ150-1.2-6.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED JEWEL MICRO LABRET		LOT NO 00972		6
4	Z-GPINMLBJ158-1.2-6.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED JEWEL MICRO LABRET		LOT NO 00972		6
5	Z-GPINMLBJ174-1.2-8.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED JEWEL MICRO LABRET		LOT NO 00972		2
6	Z-GPINMLB517PL-1.2-6.0-2.5-WOP	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED MICRO LABRET WITH PRONG SET ROUND SYNTHETIC OPAL ATTACHME		LOT NO 00972		2
7	Z-GPMBN-1.2-10.0-3	STERILE GOLD PVD COATED SURGICAL STEEL CURVED MICRO BARBELL		LOT NO 00972		4
8	Z-GPMBNSC102-1.2-8.0-3-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWEL FANCY CURVED MICRO BARBELL		LOT NO 00972		1
9	Z-GPMBNSC62-1.2-8.0-3-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWEL FANCY CURVED MICRO BARBELL - TRIPLE GEM CURVE		LOT NO 00972		2
10	Z-GPMLB-1.2-7.0-2.5	STERILE GOLD PVD COATED SURGICAL STEEL MICRO LABRET		LOT NO 00972		4
11	Z-GPPFSJ32-0.46-4.5-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWEL THREADLESS ATTACHMENT		LOT NO 00972		2
12	Z-GPSCSH005-1.2-8.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL PRONG SET CRYSTAL		LOT NO 00972		1
13	Z-GPSCSH023-1.2-8.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL CRYSTALINE JEWEL MULTI PURPOSE CLICKER		LOT NO 00972		1
14	Z-INMLBJ217-1.2-7.0-CR	STERILE SURGICAL STEEL INTERNALLY THREADED JEWEL MICRO LABRET		LOT NO 00972		2
15	Z-MBLSSC446-1.2-6.0-3-CR	STERILE SURGICAL STEEL JEWEL TRAGUS MICRO BARBELL		LOT NO 00972		2
16	Z-MBNBJ-1.2-10.0-3-CR	STERILE SURGICAL STEEL DOUBLE JEWEL CURVED MICRO BARBELL		LOT NO 00972		2
17	Z-MBNSC24-1.2-8.0-3-CR	STERILE SURGICAL STEEL JEWEL FANCY CURVED MICRO BARBELL - LIGHTNING		LOT NO 00972		2
18	Z-ZGTMBN-1.2-7.0-3	STERILE ZIRCON GOLD PVD COATED TITANIUM CURVED MICRO BARBELL		LOT NO 00972		2
19	Z-ZGTMBN-1.2-8.0-3	STERILE ZIRCON GOLD PVD COATED TITANIUM CURVED MICRO BARBELL		LOT NO 00972		4
20	Z-GPINMLBJ109-1.0-6.0-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED JEWEL MICRO LABRET		LOT NO 00972		2
22	Z-GPPFSJ33-0.46-4.5-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWEL THREADLESS ATTACHMENT - CRESCENT 3 GEMS		LOT NO 00972		4
5	Z-JB-1.5-3-CR	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL BALL	208935	LOT NO 00972		5
6	Z-JMB-CN-1.2-2.5-CR	STERILE SURGICAL STEEL VALUE JEWEL MICRO BALL		LOT NO 00972		12
7	Z-TMLB-PIN-1.2-7.0-4B-HP	STERILE TITANIUM MICRO LABRET PIN		LOT NO 00972		15
8	Z-TMLBJ-1.2-8.0-2.5-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		20
12	Z-TBNBJ-CN-1.6-10.0-5B-HP-CR	STERILE TITANIUM DOUBLE VALUE JEWEL NAVEL BANANA		LOT NO 00972		20
2	Z-BNDJ-1.5-14.0-5B-CR	STERILE SURGICAL STEEL DOUBLE PREMIUM CRYSTALS JEWEL NAVEL BANANA	208937	LOT NO 00972		72
3	Z-MBL-1.2-10.0-3	STERILE SURGICAL STEEL MICRO BARBELL		LOT NO 00972		150
4	Z-XMBN-1.1-12.10.0-3-CL	STERILE BIOFLEX MICRO CURVED BARBELL WITH STEEL BALLS		LOT NO 00972		54
1	Z-BNDJ-1.5-14.0-5B-CR	STERILE SURGICAL STEEL DOUBLE PREMIUM CRYSTALS JEWEL NAVEL BANANA	208938	LOT NO 00972		28
2	Z-XMBN-1.1-12.10.0-3-CL	STERILE BIOFLEX MICRO CURVED BARBELL WITH STEEL BALLS		LOT NO 00972		46
1	Z-GPINMLB-1.0-6.0-3	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED MICRO LABRET	208963	LOT NO 00972		1
2	Z-INMLB-1.2-8.0-2	STERILE SURGICAL STEEL INTERNALLY THREADED MICRO LABRET		LOT NO 00972		1
3	Z-XBN-1.1-6.14.0-4B-CL	STERILE BIOFLEX NAVEL BANANA WITH STEEL BALLS		LOT NO 00972		1
4	Z-XLB-1.1-5.6.0-4B-CL	STERILE BIOFLEX LABRET WITH STEEL BALL		LOT NO 00972		1
2	Z-BLC-1.6-16.0-4X4	STERILE SURGICAL STEEL BARBELL WITH CONES	208964	LOT NO 00972		1
35	Z-SCRC020-1.2-6.0	SURGICAL STEEL ROOK CLICKER - PLAIN SMOOTH	209004	LOT NO 00972		3
12	Z-GPINMLBJ-1.2-8.0-2-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWEL MICRO LABRET	209016	LOT NO 00972		40
13	Z-GPINMLBJ-1.2-8.0-3-CR	STERILE GOLD PVD COATED SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		40
14	Z-GPMCB-1.2-10.0-2.5	STERILE GOLD PVD COATED SURGICAL STEEL MICRO CIRCULAR BARBELL		LOT NO 00972		10
15	Z-GPMCB-1.2-10.0-3	STERILE GOLD PVD COATED SURGICAL STEEL MICRO CIRCULAR BARBELL		LOT NO 00972		5
16	Z-GPMCB-1.2-8.0-3	STERILE GOLD PVD COATED SURGICAL STEEL MICRO CIRCULAR BARBELL		LOT NO 00972		20
17	Z-INMLB-1.2-10.0-3	STERILE SURGICAL STEEL INTERNALLY THREADED MICRO LABRET		LOT NO 00972		20
18	Z-INMLB-1.2-8.0-3	STERILE SURGICAL STEEL INTERNALLY THREADED MICRO LABRET		LOT NO 00972		20
19	Z-INMLBJ-1.0-8.0-2-CR	STERILE SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		15
20	Z-INMLBJ-1.2-10.0-2-CR	STERILE SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		10
21	Z-INMLBJ-1.2-10.0-3-CR	STERILE SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		10
22	Z-INMLBJ-1.2-8.0-2-CR	STERILE SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		10
23	Z-INMLBJ-1.2-8.0-3-CR	STERILE SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		20
24	Z-MBN-1.2-10.0-2.5	STERILE SURGICAL STEEL CURVED MICRO BARBELL		LOT NO 00972		20
25	Z-MBN-1.2-8.0-2.5	STERILE SURGICAL STEEL CURVED MICRO BARBELL		LOT NO 00972		20
3	Z-JNO-1.0-6.5-2.35-CR	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL CURVED NOSE STUD	209035	LOT NO 00972		300
4	Z-MLB-1.2-10.0-3	STERILE SURGICAL STEEL MICRO LABRET		LOT NO 00972		100
5	Z-MLB-1.2-8.0-3	STERILE SURGICAL STEEL MICRO LABRET		LOT NO 00972		100
6	Z-TMLBJ-1.2-8.0-3-HP-CR	STERILE TITANIUM MICRO LABRET WITH JEWEL BALL		LOT NO 00972		100
1	Z-MBL-1.2-38.0-4	STERILE SURGICAL STEEL MICRO BARBELL	209097	LOT NO 00972		2
1	Z-MLB-P-1.0-10.0-3-CR	STERILE SURGICAL STEEL JEWEL MICRO LABRET	209098	LOT NO 00972		19
1	Z-MLB-P-1.0-10.0-3-CR	STERILE SURGICAL STEEL JEWEL MICRO LABRET	209099	LOT NO 00972		39
1	Z-GPPFSJ42-0.46-4.5-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWEL THREADLESS ATTACHMENT	207598	LOT NO 00972		3
3	Z-BL-1.6-29.0-5	STERILE SURGICAL STEEL BARBELL	209123	LOT NO 00972		30
7	Z-MBL-1.2-10.0-3	STERILE SURGICAL STEEL MICRO BARBELL		LOT NO 00972		200
11	Z-B-1.6-5	STERILE SURGICAL STEEL BALL	209126	LOT NO 00972		6
12	Z-BKTM-LB-1.2-9.0-3	STERILE BLACK PVD COATED TITANIUM MICRO LABRET		LOT NO 00972		4
13	Z-JMB-1.2-2.5-AB	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL MICRO BALL		LOT NO 00972		10
14	Z-JMB-1.2-2.5-CR	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL MICRO BALL		LOT NO 00972		10
15	Z-JMB-1.2-3-AB	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL MICRO BALL		LOT NO 00972		10
16	Z-JMB-1.2-3-AB	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL MICRO BALL		LOT NO 00972		10
17	Z-JMB-1.2-3-AD	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL MICRO BALL		LOT NO 00972		9
18	Z-JMB-1.2-3-BKD	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL MICRO BALL		LOT NO 00972		6
19	Z-TBL-PIN-1.6-18.0-HP	STERILE TITANIUM BARBELL PIN		LOT NO 00972		5
20	Z-TBL-PIN-1.6-19.0-HP	STERILE TITANIUM BARBELL PIN		LOT NO 00972		2
21	Z-TBNBJ-1.6-10.0-5B-HP-AQ	STERILE TITANIUM DOUBLE JEWEL NAVEL BANANA		LOT NO 00972		6
22	Z-TBNBJ-1.6-10.0-5B-HP-RO	STERILE TITANIUM DOUBLE JEWEL NAVEL BANANA		LOT NO 00972		6
23	Z-TINMLB-PIN-1.2-6.0-4B-HP	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET		LOT NO 00972		5
24	Z-TMLB-PIN-1.2-8.0-4B-HP	STERILE TITANIUM MICRO LABRET PIN		LOT NO 00972		10
25	Z-TMLBJ-1.2-8.0-2.5-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWEL MICRO LABRET		LOT NO 00972		6
2	Z-GPBCSH-1.2-12.0	STERILE GOLD PVD COATED SURGICAL STEEL HINGED SEGMENT RING	209344	LOT NO 00972		15
3	Z-TBCSH-0.8-7.0-HP	STERILE TITANIUM HINGED SEGMENT RING		LOT NO 00972		20
4	Z-TBCSH-0.8-8.0-HP	STERILE TITANIUM HINGED SEGMENT RING		LOT NO 00972		20
5	Z-TBCSH-1.0-10.0-HP	STERILE TITANIUM HINGED SEGMENT RING		LOT NO 00972		20
6	Z-TBCSH-1.0-11.0-HP	STERILE TITANIUM HINGED SEGMENT RING		LOT NO 00972		20
7	Z-TBCSH-1.0-8.0-HP	STERILE TITANIUM HINGED SEGMENT RING		LOT NO 00972		30
8	Z-TBCSH-1.0-9.0-HP	STERILE TITANIUM HINGED SEGMENT RING		LOT NO 00972		30
9	Z-TBN-PIN-1.6-10.0-HP	STERILE TITANIUM CURVED BARBELL PIN		LOT NO 00972		30
10	Z-TMLB-PIN-1.2-11.0-4B-HP	STERILE TITANIUM MICRO LABRET PIN		LOT NO 00972		20
11	Z-TMLB-PIN-1.2-9.0-4B-HP	STERILE TITANIUM MICRO LABRET PIN		LOT NO 00972		30
3	Z-MLB-P-1.2-6.0-3-CR	STERILE SURGICAL STEEL JEWEL MICRO LABRET	209137	LOT NO 00972		10
2	Z-XINLBUD-1.0-7.0-2-CL-CR	STERILE BIOFLEX INTERNAL MICRO LABRET WITH JEWEL DISC	209141	LOT NO 00972		5
1	Z-MBLJ11-1.2-8.0-3-CR	STERILE SURGICAL STEEL ROUND PRONG SET JEWEL TRAGUS MICRO BARBELL	209266	LOT NO 00972		160
3	Z-TBNBJ-1.6-12.0-5B-HP-CR	STERILE TITANIUM DOUBLE JEWEL NAVEL BANANA		LOT NO 00972		100
4	Z-TJNO-1.0-6.5-2.35-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWEL CURVED NOSE STUD		LOT NO 00972		100
5	Z-TMLB-1.2-10.0-3-HP	STERILE TITANIUM MICRO BARBELL		LOT NO 00972		100
3	Z-TBCSH-1.2-8.0-HP	STERILE TITANIUM HINGED SEGMENT RING	209327	LOT NO 00972		10
4	Z-TINML-1.6-14.0-4-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL		LOT NO 00972		8
5	Z-TINBL-1.6-34.0-4-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL		LOT NO 00972		7
6	Z-TINMBN-1.2-10.0-3-HP	STERILE TITANIUM INTERNALLY THREADED CURVED MICRO BARBELL		LOT NO 00972		10
7	Z-TINMBN-1.2-8.0-2-HP	STERILE TITANIUM INTERNALLY THREADED CURVED MICRO BARBELL		LOT NO 00972		10
8	Z-TMCB-1.2-8.0-2-HP	STERILE TITANIUM MICRO CIRCULAR BARBELL		LOT NO 00972		8
13	Z-LJNO-1.0-6.5-1.8-CR	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWEL LARGE LEFT CURVE NOSE STUD	209343	LOT NO 00972		86
14	Z-MCB-1.2-10.0-3	STERILE SURGICAL STEEL MICRO CIRCULAR BARBELL		LOT NO 00972		81
15	Z-MCB-1.2-8.0-3	STERILE SURGICAL STEEL MICRO CIRCULAR BARBELL		LOT NO 00972		16
16	Z-MLB-P-1.2-10.0-2.5-CR	STERILE SURGICAL STEEL JEWEL MICRO LABRET		LOT NO 00972		20
17	Z-MLB-P-1.2-8.0-2.5-CR	STERILE SURGICAL STEEL JEWEL MICRO LABRET		LOT NO 00972		27



26	THD-Z-TINBL-1.6-8.0-4	STERILE TITANIUM INTERNALLY THREADED BARBELL	2535447	LOT NO 00972		400
7	THD-Z-INLSCS1-1.6-5	STERILE SURGICAL STEEL ATTACHMENT FOR 1.6MM INTERNALLY THREADED PINS	2535974	LOT NO 00972		300
10	THD-Z-INSCP413-1.2	STERILE SURGICAL STEEL MICRO ATTACHMENT FOR 1.2MM INTERNALLY THREADED PINS		LOT NO 00972		300
3	TDI-Z-TJNO-OP-1.0-6.5-2.35-CR	STERILIZE CZECH CRYSTAL CURVED NOSE STUD FOR TDI	2538088	LOT NO 00972		300
12	TDI-Z10G18TMLJO-1.2-8.0-3-CR	STERILE GOLD PVD 18K COATED TITANIUM MICRO LABRET WITH OPTIMA CRYSTAL BALL		LOT NO 00972	SET/30	300
13	TDI-Z10G18TNDJO-1.5-12.0-CR	STERILE GOLD PVD 18K COATED TITANIUM DOUBLE OPTIMA CRYSTAL JEWEL NAVEL BANANA	2538841	LOT NO 00972	SET/20	200
1	Z-GPSCSS56	STERILE GOLD PVD COATED SURGICAL STEEL EAR STUDS PAIR		LOT NO 00972		10
1	ZIDBLSCC2-1.6-36.0-5	STERILE SURGICAL STEEL INDUSTRIAL BARBELL - ARROW	2539090	LOT NO 00972		20
1	PT-Z-GPTBL-1.6-16.0-5	STERILE GOLD PVD COATED TITANIUM BARBELL	2539357	LOT NO 00972		100
2	PT-Z-GPTBNJ-1.6-10.0-5/8-CR	STERILE GOLD PVD COATED TITANIUM PREMIUM CRYSTAL JEWEL NAVEL BANANA		LOT NO 00972		100
3	PT-Z-GPTBNJ-1.6-12.0-5/8-CR	STERILE GOLD PVD COATED TITANIUM PREMIUM CRYSTAL JEWEL NAVEL BANANA		LOT NO 00972		100
24	PT-Z10-TBNJ-1.6-12.0-5/8-CR	STERILE TITANIUM JEWEL NAVEL BANANA FOR PTIT ROC		LOT NO 00972	SET/50	500
1	TDI-Z-TBNJ-OP-1.6-10.0-5/8-CR	STERILE TITANIUM JEWEL NAVEL BANANA FOR TDI	2539505	LOT NO 00972		300
5	PT-Z-TBNJ-1.6-10.0-5/8-CR	STERILE TITANIUM JEWEL NAVEL BANANA FOR PTIT ROC	2539785	LOT NO 00972		100
8	PT-Z-GPTMLBJ-1.2-8.0-3-CR	STERILE GOLD PVD COATED TITANIUM PREMIUM CRYSTAL JEWEL MICRO LABRET	2538406	LOT NO 00972		50
4	TDI-Z-TBNJ-OP-1.6-12.0-5/8-CR	STERILE TITANIUM JEWEL NAVEL BANANA FOR TDI	2540575	LOT NO 00972		300
1	ZSNEEDLES-0.8-48.0	STERILIZED STAINLESS STEEL NEEDLE	2540658	LOT NO 00972		850
2	ZSNEEDLES-1.0-48.0	STERILIZED STAINLESS STEEL NEEDLE		LOT NO 00972		1300
3	Z-TIPFNUB-1.2-8.0-3-HP-CR	STERILE TITANIUM JEWEL THREADED LABRET		LOT NO 00972		2
1	TDI-Z-MLBJ-OP-1.2-8.0-3-CR	STERILE SURGICAL STEEL CZECH CRYSTAL JEWEL MICRO LABRET FOR TDI	2541028	LOT NO 00972		300
2	TDI-Z-TBNDJO-OP-1.6-10.0-5/8-CR	STERILE TITANIUM JEWEL NAVEL BANANA FOR TDI		LOT NO 00972		300
4	PT-Z-GPTMLB-1.2-8.0-3	STERILE GOLD PVD COATED TITANIUM MICRO BARBELL	2541186	LOT NO 00972		450
5	PT-Z-GPTMCB-1.2-8.0-3	STERILE GOLD PVD COATED TITANIUM MICRO CIRCULAR BARBELL		LOT NO 00972		100
9	PT-Z-GPTNO-0.8-6.5-2	STERILE GOLD PVD COATED TITANIUM CURVED BALL NOSE STUD		LOT NO 00972		100
10	PT-Z-TBCSH-0.8-8.0-HP	STERILE TITANIUM HINGED SEGMENT RING FOR PTIT ROC		LOT NO 00972		100
11	PT-Z-TBCSH-1.6-12.0-HP	STERILE TITANIUM HINGED SEGMENT RING FOR PTIT ROC		LOT NO 00972		30
12	PT-Z-TBN-1.6-12.0-4	STERILE TITANIUM TITANIUM CURVED BARBELL FOR PTIT ROC		LOT NO 00972		200
13	PT-Z-TBN-1.6-14.0-4	STERILE TITANIUM TITANIUM CURVED BARBELL FOR PTIT ROC		LOT NO 00972		50
14	PT-Z-TBNDJ-1.6-12.0-5/8-CR	STERILE TITANIUM DOUBLE JEWEL NAVEL BANANA FOR PTIT ROC		LOT NO 00972		100
15	PT-Z-TINBL-1.6-16.0-5-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL		LOT NO 00972		50
16	PT-Z-TINBN-1.2-10.0-3-HP	STERILE TITANIUM INTERNALLY THREADED CURVED MICRO BARBELL FOR PTIT ROC		LOT NO 00972		50
17	PT-Z-TJNO-0.8-6.5-2.35-CR	STERILE TITANIUM JEWEL CURVED NOSE STUD FOR PTIT ROC		LOT NO 00972		100
18	PT-Z-TMBN-1.2-14.0-3	STERILE TITANIUM MICRO CURVED BARBELL FOR PTIT ROC		LOT NO 00972		50
19	PT-Z-TMLB-1.2-7.0-3	STERILE TITANIUM MICRO LABRET FOR PTIT ROC		LOT NO 00972		150
20	PT-Z-TMLB-1.1-12.10-0-3/4B-CL	STERILE BIOFLEX MICRO LABRET WITH STEEL BALL FOR PTIT ROC		LOT NO 00972		50
21	PT-Z-TMLB-1.1-12.10-0-3/4B-CL	STERILE BIOFLEX MICRO LABRET WITH STEEL BALL FOR PTIT ROC		LOT NO 00972		100
22	PT-Z-TMLB-1.1-12.14.0-3/4B-CL	STERILE BIOFLEX MICRO LABRET WITH STEEL BALL FOR PTIT ROC		LOT NO 00972		50
23	PT-Z10-TBNDJ-1.6-12.0-5/8-CR	STERILE TITANIUM DOUBLE JEWEL NAVEL BANANA FOR PTIT ROC		LOT NO 00972	SET/50	500
24	PT-Z10-TMLB-1.2-10.0-3	STERILE TITANIUM MICRO LABRET FOR PTIT ROC		LOT NO 00972	SET/120	1200
25	PT-Z10-TMLB-1.2-8.0-3	STERILE TITANIUM MICRO LABRET FOR PTIT ROC		LOT NO 00972	SET/120	1250
26	PT-Z10-TMLBJ-1.2-10.0-2.5-CR	STERILE TITANIUM MICRO LABRET WITH JEWEL BALL FOR PTIT ROC		LOT NO 00972	SET/15	150
27	PT-Z10-TMLBJ-1.2-8.0-2.5-CR	STERILE TITANIUM MICRO LABRET WITH JEWEL BALL FOR PTIT ROC		LOT NO 00972	SET/15	150
28	PT-Z10-TMLBJ-1.2-8.0-3-CR	STERILE TITANIUM MICRO LABRET WITH JEWEL BALL FOR PTIT ROC		LOT NO 00972	SET/20	200
1	TDI-Z-TMLB-1.2-10.0-3	STERILE TITANIUM MICRO LABRET FOR TDI	2541351	LOT NO 00972		200
1	TDI-Z-TMLB-1.2-8.0-3	STERILE TITANIUM MICRO LABRET FOR TDI	2541517	LOT NO 00972		300
3	TDI-Z10TMCB-1.2-8.0-3	STERILE TITANIUM MICRO CIRCULAR BARBELL		LOT NO 00972	SET/50	500
4	TDI-Z-MBL-1.2-10.0-3	STERILE SURGICAL STEEL MICRO BARBELL FOR TDI		LOT NO 00972		100
1	TDI-Z-TBN-1.6-10.0-5/8	STERILE TITANIUM NAVEL BANANA FOR TDI	2541853	LOT NO 00972		150
1	TDI-Z-TBL-1.6-14.0-5	STERILE TITANIUM BARBELL FOR TDI	2542044	LOT NO 00972		300
2	ZSNEEDLES-1.6-48.0	STERILIZED STAINLESS STEEL NEEDLE	2542238	LOT NO 00972		3050
1	THD-Z-BKBL-1.6-12.0-4	STERILE BLACK PVD COATED SURGICAL STEEL BARBELL FOR THENDRO		LOT NO 00972		100
3	THD-Z-BKBN-1.2-8.0-3	STERILE BLACK PVD COATED SURGICAL STEEL CURVED MICRO BARBELL FOR THENDRO		LOT NO 00972		150
6	THD-Z-GPBL-1.6-16.0-5	STERILE GOLD PVD COATED SURGICAL STEEL BARBELL FOR THENDRO		LOT NO 00972		100
12	THD-Z-MLB-1.2-6.0-3	STERILIZED SURGICAL STEEL MICRO LABRET FOR THENDRO		LOT NO 00972		100
13	THD-Z-SBNDJ-1.6-10.0-4/8-CR	STERILE SURGICAL STEEL DOUBLE JEWEL MINI NAVEL BANANA FOR THENDRO		LOT NO 00972		250
16	THD-Z-TINBL-1.6-38.0-5	STERILE TITANIUM INTERNALLY THREADED BARBELL		LOT NO 00972		60
24	THD-Z-TINMC-1.2-3-HP	TITANIUM MICRO CONE FOR 1.2MM INTERNALLY THREADED PINS		LOT NO 00972		300
27	THD-Z-TINMLB-1.2-6.0-3	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET		LOT NO 00972		300
29	THD-Z-TLB-1.6-10.0-4	STERILE TITANIUM LABRET FOR THENDRO		LOT NO 00972		200
30	THD-Z-TLB-1.6-8.0-4	STERILE TITANIUM LABRET FOR THENDRO		LOT NO 00972		100
31	THD-Z-TMLB-1.2-6.0-3	STERILE TITANIUM MICRO BARBELL FOR THENDRO		LOT NO 00972		100
32	THD-Z-TMBN-1.2-6.0-3	STERILE TITANIUM MICRO CURVED BARBELL FOR THENDRO		LOT NO 00972		100
11	PT-Z-INP-1.2-50.0	STERILE SURGICAL STEEL INTERNALLY THREADED INSERTION PINS FOR PTIT ROC	2542423	LOT NO 00972		50
12	PT-Z-TBCSH-1.2-8.0-HP	STERILE TITANIUM HINGED SEGMENT RING FOR PTIT ROC		LOT NO 00972		100
15	PT-Z-TBNJ-1.6-12.0-5/8-CR	STERILE TITANIUM JEWEL NAVEL BANANA FOR PTIT ROC		LOT NO 00972		100
16	PT-Z-TINMLB-1.2-8.0-3-HP	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET FOR PTIT ROC		LOT NO 00972		200
17	PT-Z-TJNO-0.8-6.5-2.35-CR	STERILE TITANIUM JEWEL CURVED NOSE STUD FOR PTIT ROC		LOT NO 00972		400
19	PT-Z-TJNO-0.8-6.5-2.35-CR	STERILE TITANIUM LARGE CURVE JEWEL NOSE STUD FOR PTIT ROC		LOT NO 00972		300
23	PT-Z-TMLB-1.2-12.0-3	STERILE TITANIUM MICRO LABRET FOR PTIT ROC		LOT NO 00972		100
1	Z-TBCSH-0.8-8.0-HP	STERILE TITANIUM HINGED SEGMENT RING	2542543	LOT NO 00972		20
2	Z-TBCSH-0.8-9.0-HP	STERILE TITANIUM HINGED SEGMENT RING		LOT NO 00972		20
4	SV-STERILIZE	SERVICE STERILE	2538029	LOT NO 00972		20
1	SV-STERILIZE	SERVICE STERILE	2540165	LOT NO 00972		30
1	SV-STERILIZE	SERVICE STERILE	2533138	LOT NO 00972		1
3	SV-STERILIZE	SERVICE STERILE	2540091	LOT NO 00972		20
2	SV-STERILIZE	SERVICE STERILE	2540931	LOT NO 00972		10
2	SV-STERILIZE	SERVICE STERILE	2542945	LOT NO 00972		4
						21913



Sterilization:
Temperature: 44.8℃
Humidity: 20.1RH
Pressures: -1.7KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 11:1

Add ED gas
Temperature: 44.8℃
Humidity: 20.2RH
Pressures: -4.0KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 11:0

Evacuation
Temperature: 44.3℃
Humidity: 20.2RH
Pressures: -0.0KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 10:50

Keeping Temp
Temperature: 43.4℃
Humidity: 20.1RH
Pressures: -1.6KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 9:50

Keep pressure qualified
Temperature: 42.3℃
Humidity: 20.1RH
Pressures: -56.8KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 9:48

Evacuation
Temperature: 41.8℃
Humidity: 20.1RH
Pressures: -3.0KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 9:37

Heating-up:
Temperature: 23.6℃
Humidity: 20.1RH
Pressures: -3.0KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 9:0

Setting:
Temperature: 45.0℃
Humidity: 55.0RH
Pressures: -60.0KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 9:0

12/5/25

Lot 972

END

Air replacement:
Temperature: 45.0℃
Humidity: 20.2RH
Pressures: -1.6KPa
Sterilization time: 0Hour
Air replace rounds: 1
round time: 0min(s)
2025-5-12 15:29

Air replacement:
Temperature: 45.9℃
Humidity: 20.1RH
Pressures: -1.6KPa
Sterilization time: 0Hour
Air replace rounds: 2
round time: 0min(s)
2025-5-12 15:9

Air replacement:
Temperature: 45.1℃
Humidity: 20.1RH
Pressures: -1.6KPa
Sterilization time: 0Hour
Air replace rounds: 3
round time: 0min(s)
2025-5-12 14:48

Air replacement:
Temperature: 46.0℃
Humidity: 20.2RH
Pressures: -1.6KPa
Sterilization time: 0Hour
Air replace rounds: 4
round time: 0min(s)
2025-5-12 14:20

Air replacement:
Temperature: 45.5℃
Humidity: 20.1RH
Pressures: -1.6KPa
Sterilization time: 0Hour
Air replace rounds: 5
round time: 0min(s)
2025-5-12 14:7

Sterilization:
Temperature: 45.4℃
Humidity: 20.1RH
Pressures: -0.0KPa
Sterilization time: 0Hour
Air replace rounds: 5
round time: 1min(s)
2025-5-12 14:1

Sterilization:
Temperature: 45.3℃
Humidity: 20.1RH
Pressures: -0.2KPa
Sterilization time: 3Hour
Air replace rounds: 5
round time: 15min(s)
2025-5-12 11:3

