

## Report Test

**M/C :** Henan Sanqiang Medical Equipment Liability Co. Ltd

**Model :** SQ-H220


**Serial No :** HNSQ2020050063

**LOT No :** 00988

**Date :** 05 June 2025

### Determination of Sterilization Operation Specification

Item Detail	Specification	Setting	Result
Chamber Setting Temperature	40.0-55.5 °C	45.0 °C	
Vacuum pressure	40-60 kPa	60 kPa	
Vacuum Time	50-60 min	60 min	
<b>Humidification</b>			
Humidity	50.0-90.0%	20.2 %	
<b>Sterilization</b>			
Sterilization Temperature	40.0-55.5 °C	43.5 °C	
<b>Gas Injection</b>			
Final pressure	0-10 kPa	4 kPa	
Final temperature of vaporizer	40.0-55.5 °C	44.8 °C	
<b>Exposure</b>			
Temperature at start point of Exposure time period	40.0-55.5 °C	43.3 °C	
Temperature at finish point of Exposure time period	40.0-55.5 °C	46.0 °C	
Exposure Time	120-240 min	180 min	
<b>Exhaust</b>			
Exhaust Time	60-80 min	75 min	
First Vacuum pressure	0-2 kPa	1.5 kPa	
Second Vacuum pressure	0-2 kPa	1.5 kPa	
Third Vacuum pressure	0-2 kPa	1.5 kPa	
Fourth Vacuum pressure	0-2 kPa	1.5 kPa	
Fifth Vacuum pressure	0-2 kPa	1.5 kPa	
<b>Storage</b>			
Storage Time	≥ 48 hour	≥ 48 hour	
<b>Chemical Test</b>			
EO sterilization temperature 40-55°C, Time Sterilization by EO 3hours	Blue	Red	Pass
<b>Biological Test</b>			
EZ Test Biological Indicator	Green	Green	Pass

  
 (.....)  
 Mr. Jutha Liam mukda

			LOT	Color	BOX	TOTAL / PCS
		210113				
2	Z-MLB-J-1 0-8 0-3-AB	STERILE SURGICAL STEEL JEWELED MICRO LABRET	LOT NO 00988			3
		210399				
1	Z-BNDJ-1 6-11 2 5/8-CR	STERILE SURGICAL STEEL DOUBLE PREMIUM CRYSTALS JEWELED NAVEL BANANA	LOT NO 00988			4
2	Z-BNDJ-1 6-12 0 4/8-CR	STERILE SURGICAL STEEL DOUBLE PREMIUM CRYSTALS JEWELED NAVEL BANANA	LOT NO 00988			5
3	Z-MCB-1 2-9 0-3	STERILE SURGICAL STEEL MICRO CIRCULAR BARBELL	LOT NO 00988			4
4	Z-MLB-J-1 2-9 0-2 5-CR	STERILE SURGICAL STEEL JEWELED MICRO LABRET	LOT NO 00988			5
		210438				
1	Z-MLB-1 2-12 0-3	STERILE SURGICAL STEEL MICRO LABRET	LOT NO 00988			100
		210486				
2	Z-JNO-0 8-6 5-2 35-CR	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWELED CURVED NOSE STUD	LOT NO 00988			100
3	Z-MLB-1 2-8 0-3	STERILE SURGICAL STEEL MICRO LABRET	LOT NO 00988			150
		210511				
4	Z-TJNO-CN-0 8-6 5-2 35-HP-CR	STERILE TITANIUM VALUE CRYSTAL JEWELED CURVED NOSE STUD	LOT NO 00988			2
		210512				
1	Z-TB-1 6-3-HP	STERILE TITANIUM BALL	LOT NO 00988			113
		2542383				
15	THD-Z-TINBL-1 6-12 0-4	STERILE TITANIUM INTERNALLY THREADED BARBELL	LOT NO 00988			200
		2542383				
33	THD-Z-14INMJS25-1 2-CR	STERILE 14K GOLD JEWELED ATTACHMENT FOR 1 2MM INTERNALLY THREADED PINS	LOT NO 00988			60
35	THD-Z-14INMPS07-1 2	STERILE 14K GOLD ATTACHMENT FOR 1 2MM INTERNALLY THREADED PINS	LOT NO 00988			40
		2545211				
24	TDI-Z10TJNO90-0 8-6 5-2 35-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWELED 90 DEGREE NOSE STUD	LOT NO 00988		SET/30	300
		2545212				
5	TDI-Z-CNO-1 0-6 5-2 5	STERILE SURGICAL STEEL CURVED BALL NOSE STUD FOR TDI	LOT NO 00988			150
		2545341				
1	TDI-Z-TINBL-1 6-36 0-5	STERILE TITANIUM INTERNALLY THREADED BARBELL FOR TDI	LOT NO 00988			50
		2545681				
2	PT-Z-GPTBL-1 6-20 0-5	STERILE GOLD PVD COATED TITANIUM BARBELL	LOT NO 00988			50
6	PT-Z-TBCSH-1 2-10 0-HP	STERILE TITANIUM HINGED SEGMENT RING FOR PTIT ROC	LOT NO 00988			100
7	PT-Z-TBL-1 6-38 0-5	STERILIZED TITANIUM BARBELL FOR PTIT ROC	LOT NO 00988			50
8	PT-Z-TINBL-1 6-16 0-5-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL	LOT NO 00988			30
9	PT-Z-TINBL-1 6-20 0-5-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL FOR PTIT ROC	LOT NO 00988			50
12	PT-Z-TJNO-0 8-6 5-2 35-CR	STERILE TITANIUM JEWELED CURVED NOSE STUD FOR PTIT ROC	LOT NO 00988			300
14	PT-Z-TIMBN-1 2-10 0-3	STERILE TITANIUM MICRO CURVED BARBELL FOR PTIT ROC	LOT NO 00988			200
16	PT-Z-TMCB-1 2-8 0-3	STERILE TITANIUM MICRO CIRCULAR BARBELL FOR PTIT ROC	LOT NO 00988			150
17	PT-Z-TMLB-1 2-10 0-3	STERILE TITANIUM MICRO LABRET FOR PTIT ROC	LOT NO 00988			500
18	PT-Z-TMLB-1 2-8 0-3	STERILE TITANIUM MICRO LABRET FOR PTIT ROC	LOT NO 00988			400
19	PT-Z-TMLBJ-1 2-8 0-2 5-CR	STERILE TITANIUM MICRO LABRET WITH JEWELED BALL FOR PTIT ROC	LOT NO 00988			150
20	PT-Z-TST-1 6-32 0-4	STERILE TITANIUM 90 DEGREE STAPLE BARBELL FOR PTIT ROC	LOT NO 00988			10
23	PT-Z-XMBN-1 1 2-10 0-3-CL	STERILE BIOFLEX MICRO CURVED BARBELL WITH STEEL BALLS FOR PTIT ROC	LOT NO 00988			25
24	PT-Z-XMLB-1 1 2-10 0-3/4B-CL	STERILE BIOFLEX MICRO LABRET WITH STEEL BALL FOR PTIT ROC	LOT NO 00988			50
25	PT-Z-XMLB-1 1 2-14 0-3/4B-CL	STERILE BIOFLEX MICRO LABRET WITH STEEL BALL FOR PTIT ROC	LOT NO 00988			30
26	PT-Z10-TBL-1 6-16 0-5	STERILIZED TITANIUM BARBELL FOR PTIT ROC	LOT NO 00988		SET/25	250
28	PT-Z10-TBNJ-1 6-10 0-5/8-CR	STERILE TITANIUM JEWELED NAVEL BANANA FOR PTIT ROC	LOT NO 00988		SET/25	250
29	PT-Z10-TBNJ-1 6-12 0-5/8-CR	STERILE TITANIUM JEWELED NAVEL BANANA FOR PTIT ROC	LOT NO 00988		SET/50	500
30	PT-Z10-TJNO-0 8-6 5-2 35-CR	STERILE TITANIUM JEWELED CURVED NOSE STUD FOR PTIT ROC	LOT NO 00988		SET/100	1000
		2545809				
5	TDI-Z-TBL-1 6-22 0-5	STERILE TITANIUM BARBELL FOR TDI	LOT NO 00988			100
25	TDI-Z-TMLB-1 2-8 0-3	STERILE TITANIUM MICRO LABRET FOR TDI	LOT NO 00988			300
31	TDI-Z-TONR-1 0-9 0-2 5	STERILE TITANIUM OPEN NOSE RING FOR TDI	LOT NO 00988			100
		2543280				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 00988			1
		2543991				
13	SV-STERILIZE	SERVICE STERILE	LOT NO 00988			2
		2544885				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 00988			20
		2545871				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 00988			5
						5909



Sterilization:  
Temperature: 44.2t  
Humidity:20.1RH  
Pressures:-2.6KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 13:23

Add EO gas  
Temperature: 43.9t  
Humidity:20.1RH  
Pressures:-3.6KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 13:22

Evacuation  
Temperature: 44.3t  
Humidity:20.1RH  
Pressures:-0.0KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 13:18

Keeping Temp  
Temperature: 43.6t  
Humidity:20.1RH  
Pressures:-1.6KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 12:18

Keep pressure qualified  
Temperature: 42.4t  
Humidity:20.1RH  
Pressures:-58.8KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 12:18

Evacuation  
Temperature: 41.8t  
Humidity:20.1RH  
Pressures:-3.0KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 12:5

Heating-up:  
Temperature: 24.0t  
Humidity:20.1RH  
Pressures:-3.0KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 11:30

Setting:  
Temperature: 45.0t  
Humidity:55.0RH  
Pressures:-00.0KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 11:20

5/6/25

END LOT 933

Air replacement:  
Temperature: 45.1t  
Humidity:20.2RH  
Pressures:-1.6KPa  
Sterilization time:0Hour  
Air replace rounds:1  
round time:0min(s)  
2025-6-5 17:50

Air replacement:  
Temperature: 46.1t  
Humidity:20.1RH  
Pressures:-1.6KPa  
Sterilization time:0Hour  
Air replace rounds:2  
round time:0min(s)  
2025-6-5 17:30

Air replacement:  
Temperature: 45.2t  
Humidity:20.2RH  
Pressures:-1.6KPa  
Sterilization time:0Hour  
Air replace rounds:3  
round time:0min(s)  
2025-6-5 17:10

Air replacement:  
Temperature: 46.3t  
Humidity:20.1RH  
Pressures:-1.0KPa  
Sterilization time:0Hour  
Air replace rounds:4  
round time:0min(s)  
2025-6-5 16:49

Air replacement:  
Temperature: 45.4t  
Humidity:20.1RH  
Pressures:-1.6KPa  
Sterilization time:0Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 16:29

Sterilization:  
Temperature: 45.4t  
Humidity:20.2RH  
Pressures:-1.2KPa  
Sterilization time:0Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 16:23

Sterilization:  
Temperature: 44.4t  
Humidity:20.1RH  
Pressures:-3.5KPa  
Sterilization time:3Hour  
Air replace rounds:5  
round time:15min(s)  
2025-6-5 13:25

ISO 9001



ISO 14001



ISO 45001



SEDEX



CSR-DIW



FM-ISO0-009 Rev.02 (05/11/2020)