


## Report Test

M/C: Henan Sanqiang Medical Equipment Liability Co. Ltd  
 Model: SQ-H220  
 Serial No: HNSQ2020050063

LOT No: 001064  
 Date: 30 September 2025

### Determination of Sterilization Operation Specification

Item Detail	Specification	Setting	Result
Chamber Setting Temperature	40.0-55.5 °C	45.0 °C	
Vacuum pressure	40-60 kPa	60 kPa	
Vacuum Time	50-60 min	60 min	
<b>Humidification</b>			
Humidity	50.0-90.0%	20.2 %	
<b>Sterilization</b>			
Sterilization Temperature	40.0-55.5 °C	43.5 °C	
<b>Gas Injection</b>			
Final pressure	0-10 kPa	4 kPa	
Final temperature of vaporizer	40.0-55.5 °C	44.8 °C	
<b>Exposure</b>			
Temperature at start point of Exposure time period	40.0-55.5 °C	43.3 °C	
Temperature at finish point of Exposure time period	40.0-55.5 °C	46.0 °C	
Exposure Time	120-240 min	180 min	
<b>Exhaust</b>			
Exhaust Time	60-80 min	75 min	
First Vacuum pressure	0-2 kPa	1.5 kPa	
Second Vacuum pressure	0-2 kPa	1.5 kPa	
Third Vacuum pressure	0-2 kPa	1.5 kPa	
Fourth Vacuum pressure	0-2 kPa	1.5 kPa	
Fifth Vacuum pressure	0-2 kPa	1.5 kPa	
<b>Storage</b>			
Storage Time	≥ 48 hour	≥ 48 hour	
<b>Chemical Test</b>			
EO sterilization temperature 40-55°C, Time Sterilization by EO 3hours	Blue	Red	Pass
<b>Biological Test</b>			
EZ Test Biological Indicator	Green	Green	Pass

  
 (.....)  
 Mr. Jutha Liammukda  
 Research and Development Engineer



			LOT	Color	BOX	TOTAL / PCS
		214427				
183	Z-GPSCSH064-1 2-8 0-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWELED MULTI PURPOSE CLICKER	LOT NO 1084			2
		214626				
11	Z-TMLBJ-1 2-10 0-3-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWELED MICRO LABRET	LOT NO 1084			50
		214630				
13	Z-XBCS-1 6-8 0-CL	STERILE BIOFLEX SEGMENT RING	LOT NO 1084			10
14	Z-XINLBJD-1 2-8 0-2-CL-CR	STERILE BIOFLEX INTERNAL MICRO LABRET WITH JEWELED DISC	LOT NO 1084			7
		214538				
1	Z-BC-1.6-10.0-4	STERILE SURGICAL STEEL BCR	LOT NO 1084			25
		214537				
1	Z-PEBSTUDS-6-WH	STERILE PAIR OF SYNTHETIC PEARL BALL EAR STUDS	LOT NO 1084			1
		214640				
1	ZUPLST-10.0-15.0	STERILE POLYMER PIERCING CORK STOPPER	LOT NO 1084			50
		214549				
4	Z-LB-1 6-10 0-3	STERILE SURGICAL STEEL LABRET	LOT NO 1084			5
5	Z-LB-1 6-8 0-3	STERILE SURGICAL STEEL LABRET	LOT NO 1084			10
6	Z-TINB-1 6-3-HP	STERILE TITANIUM BALL FOR 1.6MM INTERNALLY THREADED PINS	LOT NO 1084			30
7	Z-TINBL-1 6-14 0-3-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL	LOT NO 1084			20
8	Z-TINBL-1 6-16 0-5-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL	LOT NO 1084			20
9	Z-TINBL-1 6-19 0-5-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL	LOT NO 1084			5
10	Z-TINBL-1 6-10 0-3-HP	STERILE TITANIUM INTERNALLY THREADED CURVED BARBELL	LOT NO 1084			5
11	Z-TINBL-1 6-12 0-4-HP	STERILE TITANIUM INTERNALLY THREADED CURVED BARBELL	LOT NO 1084			10
12	Z-TINBND-1 6-10 0-5-8-HP-CR	STERILE TITANIUM INTERNALLY THREADED DOUBLE JEWELED NAVEL BANANA BALL TOP - BOTTOM PRONG SET	LOT NO 1084			20
13	Z-TINBND-1 2-8 0-3-HP	STERILE TITANIUM INTERNALLY THREADED CURVED MICRO BARBELL	LOT NO 1084			40
14	Z-TINMLB-FW-1 2-8 0-38-HP	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET PIN FOR 0.9 MM THREAD	LOT NO 1084			20
		214651				
1	Z-TINBL-1 6-12 0-3-HP	STERILE TITANIUM INTERNALLY THREADED BARBELL	LOT NO 1084			30
3	Z-TINCB-1 6-10 0-2 5-HP	STERILE TITANIUM INTERNALLY THREADED CIRCULAR BARBELL	LOT NO 1084			30
4	Z-TINJB-1 6-3-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWELED BALL FOR 1.6MM INTERNALLY THREADED PINS	LOT NO 1084			10
5	Z-TINMLB-1 2-12 0-3-HP	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET	LOT NO 1084			400
		2556756			SET140	
2	TDI-Z10G16TMLB-1 2-8 0-3	STERILE GOLD PVD 18K COATED TITANIUM MICRO LABRET	LOT NO 1084			100
		2556895				
1	PT-Z-GPTBCSH-0 8-8 0	STERILE GOLD PVD COATED TITANIUM HINGED SEGMENT RING	LOT NO 1084			300
		2557182				
7	TDI-Z-8NJ-OP-1 6-10 0-5-8-CR	STERILE SURGICAL STEEL CZECH CRYSTAL JEWELED NAVEL BANANA FOR TDI	LOT NO 1084			200
		2557682				
1	PT-Z-GPTMLB-1 2-8 0-3	STERILE GOLD PVD COATED TITANIUM MICRO LABRET	LOT NO 1084			200
		2557769				
2	TDI-Z-TNO-1 0-6 5-2	STERILE TITANIUM TITANIUM CURVED BALL NOSE STUD FOR TDI	LOT NO 1084			30
		2557858				
1	Z-LINO-1 0-8 0-2 95-CR	STERILE SURGICAL STEEL PREMIUM CRYSTAL JEWELED LARGE LEFT CURVE NOSE STUD	LOT NO 1084			10
2	Z-TINCB-1 6-8 0-4-HP	STERILE TITANIUM INTERNALLY THREADED CIRCULAR BARBELL	LOT NO 1084			1660

Sterilization:  
Temperature: 44.6°C  
Humidity: 20.1RH  
Pressures: -0.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 13:39

Add EO gas  
Temperature: 44.4°C  
Humidity: 20.2RH  
Pressures: -1.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 13:39

Evacuation  
Temperature: 44.2°C  
Humidity: 20.1RH  
Pressures: -0.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 13:34

Keeping Temp  
Temperature: 42.8°C  
Humidity: 20.1RH  
Pressures: -1.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 12:34

Keep pressure qualified  
Temperature: 41.5°C  
Humidity: 20.2RH  
Pressures: -58.7KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 12:33

Evacuation  
Temperature: 41.6°C  
Humidity: 20.2RH  
Pressures: -3.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 12:21

Heating-up:  
Temperature: 23.1°C  
Humidity: 20.1RH  
Pressures: -3.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 11:41

Setting:  
Temperature: 45.0°C  
Humidity: 55.0RH  
Pressures: -0.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 11:40

30/09/25  
END Lot 4064

Air replacement:  
Temperature: 45.4°C  
Humidity: 20.2RH  
Pressures: -1.0KPa  
Sterilization time: 0Hour  
Air replace rounds: 1  
round time: 0min(s)  
2025-9-29 18:7

Air replacement:  
Temperature: 45.3°C  
Humidity: 20.1RH  
Pressures: -1.0KPa  
Sterilization time: 0Hour  
Air replace rounds: 2  
round time: 0min(s)  
2025-9-29 17:47

Air replacement:  
Temperature: 45.2°C  
Humidity: 20.2RH  
Pressures: -1.0KPa  
Sterilization time: 0Hour  
Air replace rounds: 3  
round time: 0min(s)  
2025-9-29 17:26

Air replacement:  
Temperature: 45.7°C  
Humidity: 20.2RH  
Pressures: -1.0KPa  
Sterilization time: 0Hour  
Air replace rounds: 4  
round time: 0min(s)  
2025-9-29 17:0

Air replacement:  
Temperature: 45.6°C  
Humidity: 20.1RH  
Pressures: -1.0KPa  
Sterilization time: 0Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 16:46

Sterilization:  
Temperature: 44.3°C  
Humidity: 20.1RH  
Pressures: -0.0KPa  
Sterilization time: 0Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 16:40

Sterilization:  
Temperature: 44.6°C  
Humidity: 20.1RH  
Pressures: -0.0KPa  
Sterilization time: 3Hour  
Air replace rounds: 5  
round time: 15min(s)  
2025-9-29 13:41