


Report Test

M/C: Henan Sanqiang Medical Equipment Liability Co. Ltd
 Model: SQ-H220
 Serial No: HNSQ2020050063

LOT No: 001144
 Date: 09 February 2026

Determination of Sterilization Operation Specification

Item Detail	Specification	Setting	Result
Chamber Setting Temperature	40.0-55.5 °C	45.0 °C	
Vacuum pressure	40-60 kPa	60 kPa	
Vacuum Time	50-60 min	60 min	
Humidification			
Humidity	50.0-90.0%	20.2 %	
Sterilization			
Sterilization Temperature	40.0-55.5 °C	43.5 °C	
Gas Injection			
Final pressure	0-10 kPa	4 kPa	
Final temperature of vaporizer	40.0-55.5 °C	44.8 °C	
Exposure			
Temperature at start point of Exposure time period	40.0-55.5 °C	43.3 °C	
Temperature at finish point of Exposure time period	40.0-55.5 °C	46.0 °C	
Exposure Time	120-240 min	180 min	
Exhaust			
Exhaust Time	60-80 min	75 min	
First Vacuum pressure	0-2 kPa	1.5 kPa	
Second Vacuum pressure	0-2 kPa	1.5 kPa	
Third Vacuum pressure	0-2 kPa	1.5 kPa	
Fourth Vacuum pressure	0-2 kPa	1.5 kPa	
Fifth Vacuum pressure	0-2 kPa	1.5 kPa	
Storage			
Storage Time	≥ 48 hour	≥ 48 hour	
Chemical Test			
EO sterilization temperature 40-55°C, Time Sterilization by EO 3hours	Blue	Red	Pass
Biological Test			
EZ Test Biological Indicator	Green	Green	Pass


 (.....)
 Mr. Jutha Liamukda
 Research and Development Engineer



FM-ISO0-009 Rev.02 (05/11/2020)

			LOT	Color	BOX	TOTAL / PCS
		218241				
20	Z-BKTMCB-1.2-9.0-3	STERILE BLACK PVD COATED TITANIUM MICRO CIRCULAR BARBELL	LOT NO 001144			2
21	Z-PFINJD-1.2-8.0-2.6-CR	STERILE SURGICAL STEEL JEWELED THREADLESS DISC LABRET	LOT NO 001144			10
22	Z-TBL-PIN-1.6-14.0-HP	STERILE TITANIUM BARBELL PIN	LOT NO 001144			10
23	Z-TBL-PIN-1.6-15.0-HP	STERILE TITANIUM BARBELL PIN	LOT NO 001144			8
24	Z-TBL-PIN-1.6-19.0-HP	STERILE TITANIUM BARBELL PIN	LOT NO 001144			8
25	Z-TBNDJ-1.6-10.0-5/8-HP-AQ	STERILE TITANIUM DOUBLE JEWELED NAVEL BANANA	LOT NO 001144			5
26	Z-TBNDJ-OP-1.6-10.0-5/8-HP-CR	STERILE TITANIUM DOUBLE OPTIMA CRYSTAL JEWELED NAVEL BANANA	LOT NO 001144			8
27	Z-TJMB-OP-1.2-3-HP-CR	STERILE TITANIUM OPTIMA CRYSTAL JEWELED MICRO BALL	LOT NO 001144			8
28	Z-XTBL-PINS-1.6-12.0-CL	STERILE BIOFLEX THREADED BARBELL PIN	LOT NO 001144			4
29	Z-XTBL-PINS-1.6-14.0-CL	STERILE BIOFLEX THREADED BARBELL PIN	LOT NO 001144			4
30	Z-XTBL-PINS-1.6-16.0-CL	STERILE BIOFLEX THREADED BARBELL PIN	LOT NO 001144			2
32	Z-JBS-1.6-4-AM	STERILE SURGICAL STEEL SIDE THREADED PREMIUM CRYSTAL JEWELED BALL	LOT NO 001144			4
		218247				
1	Z-TMLB-1.2-12.0-3/4-HP	STERILE TITANIUM MICRO LABRET	LOT NO 001144			50
2	Z-TMLB-1.2-10.0-3-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWELED MICRO LABRET	LOT NO 001144			50
		218256				
1	Z-TBL-1.6-38.0-4-HP	STERILE TITANIUM BARBELL	LOT NO 001144			10
		218260				
5	ZDPU-2.5	STERILE DERMAL PUNCH	LOT NO 001144			10
8	ZDPU-2.0	STERILE DERMAL PUNCH	LOT NO 001144			20
		218264				
2	Z-GPTMGB-1.2-10.0-3	STERILE GOLD PVD COATED TITANIUM MICRO CIRCULAR BARBELL	LOT NO 001144			10
3	Z-TBL-1.6-8.0-3-HP	STERILE TITANIUM BARBELL	LOT NO 001144			15
		218266				
1	Z-SCNS53-1.6-16.0	STERILE SURGICAL STEEL NIPPLE BAR - WRENCH	LOT NO 001144			1
		218306				
1	Z-TINMLBJ-1.0-8.0-2-HP-CR	STERILE TITANIUM INTERNALLY THREADED JEWELED MICRO LABRET	LOT NO 001144			20
2	Z-TJNO-1.0-6.5-2.35-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWELED CURVED NOSE STUD	LOT NO 001144			20
		2568223				
2	Z-SCBCSH-3.0	STERILE SURGICAL STEEL HINGED SEGMENT RING	LOT NO 001144			19
		2568539				
1	Z-TINMBL-1.2-10.0-3-HP	STERILE TITANIUM INTERNALLY THREADED MICRO BARBELL	LOT NO 001144			100
		2570160				
12	PT-Z-TBL-1.6-12.0-4	STERILIZED TITANIUM BARBELL FOR PTIT ROC	LOT NO 001144			100
		2570335				
1	THD-Z-BKBC-1.2-8.0-4	STERILE BLACK PVD COATED SURGICAL STEEL BALL CLOSURE RING FOR THENDRO	LOT NO 001144			100
2	THD-Z-BKBL-1.6-16.0-5	STERILE BLACK PVD COATED SURGICAL STEEL BARBELL FOR THENDRO	LOT NO 001144			100
5	THD-Z-BKLB-1.6-8.0-4	STERILE BLACK PVD COATED SURGICAL STEEL LABRET FOR THENDRO	LOT NO 001144			100
7	THD-Z-BKLAB-1.2-12.0-3	STERILE BLACK PVD COATED SURGICAL STEEL MICRO LABRET FOR THENDRO	LOT NO 001144			100
32	THD-Z-JNO-1.0-6.5-2.35-AM	STERILE SURGICAL STEEL JEWELED CURVED NOSE STUD FOR THENDRO	LOT NO 001144			120
34	THD-Z-JNO-1.0-6.5-2.35-LR	STERILE SURGICAL STEEL JEWELED CURVED NOSE STUD FOR THENDRO	LOT NO 001144			100
41	THD-Z-TIFSTD-1.6-10.0-4	STERILE FLAT INTERNALLY THREADED TITANIUM SURFACE BARBELL WITH DISCS	LOT NO 001144			100
		2570857				
13	PT-Z10-TBNDJ-1.6-12.0-5/8-CR	STERILE TITANIUM DOUBLE JEWELED NAVEL BANANA FOR PTIT ROC	LOT NO 001144		SET/50	500
		2571386				
1	Z-GPTBNDJ-1.6-12.0-5/8-CR	STERILE GOLD PVD COATED TITANIUM DOUBLE PREMIUM CRYSTAL JEWELED NAVEL BANANA	LOT NO 001144			10
		2571428				
2	TDI-Z-TINMLB-1.2-6.0-2	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET FOR TDI	LOT NO 001144			100
3	TDI-Z-TINMLB-1.2-6.0-3	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET FOR TDI	LOT NO 001144			50
4	TDI-Z-TINMLB-1.2-8.0-2	STERILE TITANIUM INTERNALLY THREADED MICRO LABRET FOR TDI	LOT NO 001144			900
		2571201				
1	ZSNEEDLES-4.0-48.0	STERILIZED STAINLESS STEEL NEEDLE	LOT NO 001144			403
		2571683				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 001144			1
		2567015				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 001144			16
						2856

Sterilization:
Temperature: 44.7°C
Humidity: 20.18%
Pressure: -4.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:17

Add ED gas
Temperature: 44.5°C
Humidity: 20.19%
Pressure: -4.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:17

Evacuation
Temperature: 44.2°C
Humidity: 20.19%
Pressure: -0.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:13

Keeping Temp
Temperature: 42.5°C
Humidity: 20.19%
Pressure: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:13

Keep pressure qualified
Temperature: 41.6°C
Humidity: 20.19%
Pressure: -0.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:11

Evacuation
Temperature: 41.5°C
Humidity: 20.19%
Pressure: -0.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:09

heating-up:
Temperature: 21.0°C
Humidity: 20.19%
Pressure: 0.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:15

Setting:
Temperature: 45.0°C
Humidity: 20.19%
Pressure: -0.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:15

END

9/2/26
Lot 1144

Air replacement:
Temperature: 45.2°C
Humidity: 20.19%
Pressure: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:44

Air replacement:
Temperature: 45.0°C
Humidity: 20.19%
Pressure: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:24

Air replacement:
Temperature: 45.5°C
Humidity: 20.23%
Pressure: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:4

Air replacement:
Temperature: 46.3°C
Humidity: 20.19%
Pressure: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 4
round time: 15min(s)
2026-2-9 16:43

Air replacement:
Temperature: 45.3°C
Humidity: 20.19%
Pressure: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 16:23

Sterilization:
Temperature: 45.1°C
Humidity: 20.19%
Pressure: -2.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 16:18

Sterilization:
Temperature: 44.2°C
Humidity: 20.19%
Pressure: -5.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-2-9 13:10

