


Report Test

M/C: Henan Sanqiang Medical Equipment Liability Co. Ltd
 Model: SQ-H220
 Serial No: HNSQ2020050063

LOT No: 001175
 Date: 30 March 2026

Determination of Sterilization Operation Specification

Item Detail	Specification	Setting	Result
Chamber Setting Temperature	40.0-55.5 °C	45.0 °C	
Vacuum pressure	40-60 kPa	60 kPa	
Vacuum Time	50-60 min	60 min	
Humidification			
Humidity	50.0-90.0%	20.2 %	
Sterilization			
Sterilization Temperature	40.0-55.5 °C	43.5 °C	
Gas Injection			
Final pressure	0-10 kPa	4 kPa	
Final temperature of vaporizer	40.0-55.5 °C	44.8 °C	
Exposure			
Temperature at start point of Exposure time period	40.0-55.5 °C	43.3 °C	
Temperature at finish point of Exposure time period	40.0-55.5 °C	46.0 °C	
Exposure Time	120-240 min	180 min	
Exhaust			
Exhaust Time	60-80 min	75 min	
First Vacuum pressure	0-2 kPa	1.5 kPa	
Second Vacuum pressure	0-2 kPa	1.5 kPa	
Third Vacuum pressure	0-2 kPa	1.5 kPa	
Fourth Vacuum pressure	0-2 kPa	1.5 kPa	
Fifth Vacuum pressure	0-2 kPa	1.5 kPa	
Storage			
Storage Time	≥ 48 hour	≥ 48 hour	
Chemical Test			
EO sterilization temperature 40-55°C, Time Sterilization by EO 3hours	Blue	Red	Pass
Biological Test			
EZ Test Biological Indicator	Green	Green	Pass


 (.....)
 Mr. Jutha Liammukda
 Research and Development Engineer



			LOT	Color	BOX	TOTAL / PCS
		219575				
3	Z-TJNO-1.0-6.5-2.35-HP-CR	STERILE TITANIUM PREMIUM CRYSTAL JEWELED CURVED NOSE STUD	LOT NO 1176			9
		219586				
4	Z-TFFIN-PIN-1.0-8.0-2.58-HP	STERILE TITANIUM THREADLESS LABRET PIN	LOT NO 1176			5
		219663				
1	Z-GPBNJ-1.6-12.0-5.8-CR	STERILE GOLD PVD COATED SURGICAL STEEL PREMIUM CRYSTAL JEWELED NAVEL BANANA	LOT NO 1176			31
2	Z-GPMLBJ-1.2-8.0-3-CR	STERILE GOLD PVD COATED SURGICAL STEEL JEWELED MICRO LABRET	LOT NO 1176			50
3	Z-IMLBJ-1.2-8.0-2-CR	STERILE SURGICAL STEEL INTERNALLY THREADED PREMIUM CRYSTAL JEWELED MICRO LABRET	LOT NO 1176			50
4	Z-SES1-2.5-CR	STERILE SURGICAL STEEL ROUND PRONG SET JEWELED EAR STUDS PAIR	LOT NO 1176			50
		219684				
1	Z-BCSH-1.2-10.0	316L STERILE SURGICAL STEEL HINGED SEGMENT RING POLISHED SEAMLESS LOOK CLICKER	LOT NO 1176			4
2	Z-BNN-1.6-11.2-4.6	STERILE SURGICAL STEEL NAVEL BANANA	LOT NO 1176			4
3	Z-BNN-1.6-14.0-5.8	STERILE SURGICAL STEEL NAVEL BANANA	LOT NO 1176			4
4	Z-CB-1.6-10.0-3	STERILE SURGICAL STEEL CIRCULAR BARBELL	LOT NO 1176			4
5	Z-PFINJD-1.0-8.0-2.5-CR	STERILE SURGICAL STEEL JEWELED THREADLESS DISC LABRET	LOT NO 1176			9
6	Z-PFINJD-1.2-8.0-3-CR	STERILE SURGICAL STEEL JEWELED THREADLESS DISC LABRET	LOT NO 1176			9
7	Z-PFIN-PIN-1.0-6.0-2.55	STERILE SURGICAL STEEL THREADLESS LABRET PIN	LOT NO 1176			6
		219710				
1	Z-TINBJ-1.6-12.0-6.8-HP-CR	STERILE TITANIUM INTERNALLY THREADED PREMIUM CRYSTAL JEWELED NAVEL BANANA	LOT NO 1176			20
		219715				
1	Z-TMBL-1.2-8.0-3-HP	STERILE TITANIUM MICRO BARBELL	LOT NO 1176			250
		219738				
4	Z-TMBN-PIN-1.2-8.0-HP	STERILE TITANIUM MICRO CURVED BARBELL PIN	LOT NO 1176			50
5	Z-TMBN-PIN-1.2-9.0-HP	STERILE TITANIUM MICRO CURVED BARBELL PIN	LOT NO 1176			50
		219744				
4	Z-BL-1.6-29.0-5	STERILE SURGICAL STEEL BARBELL	LOT NO 1176			5
		219752				
1	Z-IMLBJD-1.0-8.0-2.5-CR	STERILE SURGICAL STEEL INTERNALLY THREADED LABRET WITH JEWELED DISC	LOT NO 1176			50
2	Z-MBLSJ1-1.2-8.0-3-CR	STERILE SURGICAL STEEL ROUND PRONG SET JEWELED TRAGUS MICRO BARBELL	LOT NO 1176			75
		2571073				
1	Z-GPBCSH-0.8-5.0	GOLD PVD COATED SURGICAL STEEL HINGED SEGMENT RING PLAIN POLISHED	LOT NO 1176			30
		257319				
3	TDI-Z-TBNJ-OP-1.6-10.0-5.8-RO	STERILE TITANIUM JEWELED NAVEL BANANA FOR TDI	LOT NO 1176			100
		2575810				
1	PT-2-GPTBL-1.6-14.0-5	STERILE GOLD PVD COATED TITANIUM BARBELL	LOT NO 1176			50
3	PT-2-GPTBL-1.6-22.0-4	STERILE GOLD PVD COATED TITANIUM BARBELL	LOT NO 1176			50
		2576681				
2	PT-2-GPTMCB-1.2-10.0-3	STERILE GOLD PVD COATED TITANIUM MICRO CIRCULAR BARBELL	LOT NO 1176			50
3	PT-2-GPTMLB-1.2-8.0-3	STERILE GOLD PVD COATED TITANIUM MICRO LABRET	LOT NO 1176			200
4	PT-2-TBL-1.6-12.0-4	STERILIZED TITANIUM BARBELL FOR PTIT ROC	LOT NO 1176			50
6	PT-2-TGN-1.6-12.0-4	STERILE TITANIUM CURVED BARBELL FOR PTIT ROC	LOT NO 1176			100
7	PT-2-TMBN-1.2-10.0-3-HP	STERILE TITANIUM INTERNALLY THREADED CURVED MICRO BARBELL FOR PTIT ROC	LOT NO 1176			50
11	PT-Z10-TMBL-1.2-8.0-3-HP	STERILE TITANIUM MICRO BARBELL FOR PTIT ROC	LOT NO 1176		SET/50	500
12	PT-Z10-TMBN-1.2-10.0-3	STERILE TITANIUM MICRO CURVED BARBELL FOR PTIT ROC	LOT NO 1176		SET/25	250
13	PT-Z10-TMLB-1.2-10.0-3	STERILE TITANIUM MICRO LABRET FOR PTIT ROC	LOT NO 1176		SET/100	1000
		2576748				
1	TDI-Z-BL-1.6-32.0-5	STERILE SURGICAL STEEL BARBELL FOR TDI	LOT NO 1176			50
2	TDI-Z-BL-1.6-38.0-5	STERILE SURGICAL STEEL BARBELL FOR TDI	LOT NO 1176			50
		2577058				
1	TDI-Z-TMBL-1.2-10.0-3-HP	STERILE TITANIUM MICRO BARBELL FOR TDI	LOT NO 1176			200
		2575000				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 1176			1
		2577193				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 1176			22
		2577374				
1	SV-STERILIZE	SERVICE STERILE	LOT NO 1176			9
						3494

Sterilization:
Temperature: 44.3°C
Humidity: 20.1RH
Pressures: -7.3kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:29

Add EO gas
Temperature: 44.1°C
Humidity: 20.2RH
Pressures: -4.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:30

Evacuation
Temperature: 44.1°C
Humidity: 20.1RH
Pressures: -8.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:34

Keeping Temp
Temperature: 42.5°C
Humidity: 20.1RH
Pressures: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:34

Keep pressure qualified
Temperature: 41.8°C
Humidity: 20.1RH
Pressures: -53.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:37

Evacuation
Temperature: 41.6°C
Humidity: 20.2RH
Pressures: -3.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:39

Heating up:
Temperature: 32.8°C
Humidity: 20.1RH
Pressures: -5.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:41

Setting:
Temperature: 45.0°C
Humidity: 20.0RH
Pressures: -8.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:43

3013126
Lot 1175

Air replacement:
Temperature: 44.3°C
Humidity: 20.1RH
Pressures: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 1
round time: 15min(s)
2026-3-30 13:47

Air replacement:
Temperature: 44.3°C
Humidity: 20.1RH
Pressures: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 2
round time: 15min(s)
2026-3-30 13:47

Air replacement:
Temperature: 45.5°C
Humidity: 20.1RH
Pressures: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 3
round time: 15min(s)
2026-3-30 13:47

Air replacement:
Temperature: 45.4°C
Humidity: 20.1RH
Pressures: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 4
round time: 15min(s)
2026-3-30 13:48

Air replacement:
Temperature: 44.5°C
Humidity: 20.2RH
Pressures: -1.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:48

Sterilization:
Temperature: 45.0°C
Humidity: 20.1RH
Pressures: -9.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:49

Sterilization:
Temperature: 44.1°C
Humidity: 20.1RH
Pressures: -0.0kPa
Sterilization time: 9hour
Air replace rounds: 5
round time: 15min(s)
2026-3-30 13:49